

# The Shoeshine



For the ultimate in self-administered shoeshine appliances, there is no finer device than the Dinklaker Shoe Attendant. The wearer of the target shoe simply steps onto the application platform, slips the Brainwave Polish Monitor (pat. pend.) over her head, then, in her mind's eye, imagines the type of shine she desires. A sophisticated artificial intelligence module collects and interprets the image, then activates the polisher, wax applicator, sock deodorizer, and buffer as necessary. If the user keeps her feet still, she will realize a shoe custom-shined to perfection every time! (Note: Dismemberment can occur if the target shoe moves even one centimeter during the deodorizing procedure, so the user's ability to remain immobile is key to a safe and pleasant outcome.)